


Laser pointer headband for testing cervical joint position sense

You can purchase lasers from <https://www.alphasport.com.au/product/laser-headlamp/>



Or

<http://www.trackerlaser.com://www.motionguidance.com>



<https://concussionlab.com/collections/all>



You can also make your own for about \$10



Purchase the lasers from ebay- search laser pointer key chain. Red lasers are the best. Only buy from Australia so that the lasers should be less than 1milliwatt which is the legal requirement. You have to have a legitimate reason for carrying a laser around. Buy the headbands from the supermarket ebay and then attach the self adhesive cord clamps from Office Works on to it. The laser then slips into the clip and can be rotated to turn and hold the laser on and then turned again to switch it off.

<http://www.officeworks.com.au/shop/officeworks/p/wrapt-up-single-cord-clamps-6-pack-crcmccp6>

